

# Product Information

## VACUETTE® Tubes Shelf Life

VACUETTE® Blood Collection Tubes	Code	Fill Volume	Shelf life
Coagulation sodium citrate 3.2%, 3,8% sandwich (SW)	9NC 3.2%	1ml	6 months
Coagulation sodium citrate 3.2%, 3,8% sandwich (SW)	9NC 3.2%, 9NC 3.8%	2 - 3.5ml	12 months
Coagulation sodium citrate 3.2%, 3,8% single-wall	9NC 3.2%, 9NC 3.8%	depending on item number*	
K2EDTA, K3EDTA spray dried	K2E, K3E	1ml	12 months
K2EDTA, K3EDTA spray dried	K2E, K3E	> 1ml	18 months
K2EDTA Sep spray dried	K2E	depending on item number*	
CTAD, sandwich	-	> 1ml	12 months
ESR sodium citrate; 13x75mm	4NC 3.2% (PET)	≥ 1ml	9 months
ESR sodium citrate; 9x120mm	4NC 3.2% (PP)	≥ 1ml	12 months
ACD, CPDA	-	> 1ml	9 months
No Additive	Z	> 1ml	18 months
Z Serum Clot Activator	Z	1ml	12 months
Z Serum Clot Activator	Z	> 1ml	18 months
Z Serum Sep Clot Activator	Z	depending on item number*	
Lithium Heparin	LH	1ml	12 months
Lithium/ Ammonium/ Sodium Heparin (Sep)	LH, LH/Gel, AH, NH	> 1ml	18 months
Sodium Fluoride / Potassium Oxalate	FX	> 1ml	18 months
EDTA Aprotinin	-	≥ 1ml	6 months
GLUCOMEDICS, liquid, sandwich	-	≥ 1ml	6 months
Homocysteine Detection Tube, sandwich	-	≥ 1ml	12 months

\* please see individual label

Printed copies will not systematically be updated. Source: SAP 19.11.2015; this file replaces all previously issued documents regarding shelf-life.  
Valid at the date of issue. Information is subject to change.

# Product Information

## VACUETTE® and MiniCollect® Shelf Life

VACUETTE® Blood Collection Accessories	Shelf life
VACUETTE® Multiple Use Drawing Needle	60 months*
VACUETTE® VISIO PLUS Needle	60 months*
VACUETTE® Luer Adapter	60 months*
Blood Collection/Infusion Sets	36 months
SAFETY Blood Collection/Infusion Sets	36 months
Blood Transfer Unit	36 months
HOLDEX® Single Use Holder	36 months
VACUETTE® QUICKSHIELD Safety Tube Holder	36 months
VACUETTE® QUICKSHIELD Complete PLUS	36 months
VACUETTE® TUBE-TOUCH	12 months
VACUETTE® Disposable Tourniquet	24 months

MiniCollect® Capillary Blood Collection System	Code	Fill Volume	Shelf life
Z Serum	Z	0.5ml, 1ml	18 months
Z Serum Sep	Z	0.5ml, 0.8ml	18 months
K2EDTA	K2E	0.25ml, 0.5ml	18 months
K3EDTA	K3E	0.25ml, 0.5ml, 1ml	18 months
Lithium Heparin	LH	0.5ml, 1ml	18 months
Lithium Heparin Sep	LH	0.8ml	18 months
Sodium Fluoride / Potassium Oxalate	FX	0.25ml, 0.5ml	18 months
Coagulation sodium citrate 3.8%	9NC 3.8%	1ml	15 months
Coagulation sodium citrate 3.2%	9NC 3.2%	0.5ml, 1ml	15 months
No Additive	Z		60 months

MiniCollect® Accessories	Shelf life
Safety Lancets, sterile	36 months

Printed copies will not systematically be updated. Source: SAP 19.11.2015; this file replaces all previously issued documents regarding shelf-life.  
Valid at the date of issue. Information is subject to change.